

Experiment File Path: C:\Users\User\Documents\Experiments\Lubomir\IFN-y_mouse_samples_Lubomir_Bedn...

Protocol File Path: C:\Users\User\Documents\Protocols\IFN-y_mouse_samples_Lubomir_Bednar.prt

Plate Number
Date
Time
Reader Type: Synergy HTX
Reader Serial Number: Unknown
Reading Type: Reader

Procedure Details

Plate Type 96 WELL PLATE (Use plate lid)
Eject plate on completion
Read IFN_read
Absorbance Endpoint
Full Plate
Wavelengths: 450
Read Speed: Normal, Delay: 100 msec, Measurements/Data Point: 8
Set Temperature Incubator off

Layout

	1	2	3	4	5	6	
A	STD1 250	STD1 250					Well ID Conc/Dil
B	STD2 125	STD2 125			SPL2:1 0,5	SPL2:1 0,5	Well ID Conc/Dil
C	STD3 62,5	STD3 62,5	SPL1:1 0,25	SPL1:1 0,25	SPL2:2 0,25	SPL2:2 0,25	Well ID Conc/Dil
D	STD4 31,25	STD4 31,25	SPL1:2 0,125	SPL1:2 0,125	SPL2:3 0,125	SPL2:3 0,125	Well ID Conc/Dil
E	STD5 15,625	STD5 15,625	SPL1:3 0,0625	SPL1:3 0,0625	SPL2:4 0,0625	SPL2:4 0,0625	Well ID Conc/Dil
F	STD6 7,8125	STD6 7,8125	SPL1:4 0,03125	SPL1:4 0,03125	SPL2:5 0,03125	SPL2:5 0,03125	Well ID Conc/Dil
G	STD7 3,9063	STD7 3,9063	SPL1:5 0,015625	SPL1:5 0,015625	SPL2:6 0,015625	SPL2:6 0,015625	Well ID Conc/Dil
H	BLK	BLK	SPL1:6 0,0078125	SPL1:6 0,0078125	SPL2:7 0,0078125	SPL2:7 0,0078125	Well ID Conc/Dil

Layout (continued)

	7	8	9	10	11	12	
A	SPL3:1 1	SPL3:1 1			STDB1 250	STDB1 250	Well ID
B	SPL3:2 0,5	SPL3:2 0,5			STDB2 125	STDB2 125	Conc/Dil Well ID
C	SPL3:3 0,25	SPL3:3 0,25	SPL4:1 0,25	SPL4:1 0,25	STDB3 62,5	STDB3 62,5	Conc/Dil Well ID
D	SPL3:4 0,125	SPL3:4 0,125	SPL4:2 0,125	SPL4:2 0,125	STDB4 31,25	STDB4 31,25	Conc/Dil Well ID
E	SPL3:5 0,0625	SPL3:5 0,0625	SPL4:3 0,0625	SPL4:3 0,0625	STDB5 15,625	STDB5 15,625	Conc/Dil Well ID
F	SPL3:6 0,03125	SPL3:6 0,03125	SPL4:4 0,03125	SPL4:4 0,03125	STDB6 7,8125	STDB6 7,8125	Conc/Dil Well ID
G	SPL3:7 0,015625	SPL3:7 0,015625	SPL4:5 0,015625	SPL4:5 0,015625	STDB7 3,9063	STDB7 3,9063	Conc/Dil Well ID
H	SPL3:8 0,0078125	SPL3:8 0,0078125	SPL4:6 0,0078125	SPL4:6 0,0078125	BLK	BLK	Conc/Dil Well ID

IFN_read:450

	1	2	3	4	5	6	7	8	9	10	11	12	
A													IFN_read:450
B													IFN_read:450
C													IFN_read:450
D													IFN_read:450
E													IFN_read:450
F													IFN_read:450
G													IFN_read:450
H													IFN_read:450

Blank IFN_read:450

	1	2	3	4	5	6	7	8	9	10	11	12	
A													Blank IFN_read:450
B													Blank IFN_read:450
C													Blank IFN_read:450
D													Blank IFN_read:450
E													Blank IFN_read:450
F													Blank IFN_read:450
G													Blank IFN_read:450
H													Blank IFN_read:450